Uri-Color Check Tiras de Urina

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ENGLISH

CLINICAL IMPORTANCE

Urine testing provides the clinician with information on possible pathologies of the urinary tract, as well as about some extra-reputable diseases. Due to its simplicity, low cost and the ease of obtaining the sample for analysis, it is a routine examination, already used for many years.

The **Uri-Color Check** (Urine Strip) of **WAMA Diagnóstica** determines semiquantitatively the presence of urobilinogen, glucose, ketones, bilirubin, protein, nitrite, pH, blood, specific gravity and leukocytes in urine samples, with a high degree of accuracy.

PRINCIPLE OF METHOD

Urobilinogen: Based on the diazotization reaction of 4-Methoxybenzene diazonium with Urobilinogen, in the presence of strong acid, the color scale of the parameter varies from light pink to dark pink, according to the concentration in the urine sample.

Glucose: This is a sequential double enzymatic reaction. In the first stage, glucose oxidase catalyzes the reaction between glucose present in the urine sample and oxygen from ambient air, producing gluconic acid and hydrogen peroxide. In the second stage, peroxidase catalyzes the reaction between the previously formed hydrogen peroxide and the potassium iodide chromogen, giving rise to a blue colored compound for negative samples and varies from green to brown for positive samples.

Ketones: Acetoacetic Acid, if present in the urine sample, reacts with sodium nitropusside in alkaline medium, giving rise to a violet color proportional to the concentration of the acid in positive samples. In its absence the parameter assumes a beige color, which characterizes negative samples.

Bilirubin: Bilirubin, if present in the urine sample, reacts in acidic medium with an organic diazonian salt, forming a color ranging from yellowish brown in negative samples to violet, which indicates the highest concentration in positive samples.

Protein: Protein, if present in the urine sample, promotes an ionic exchange with the blue indicator of tetrabromophenol, without causing major changes in the pH of the medium. The color scale varies from light yellow, in negative samples, to intense green, in positive samples.

Nitrite: Nitrite is a product of nitrate reduction, a normal constituent of urine, carried out by certain bacteria. In the samples, nitrite reacts with an organic acid, forming a diazonium salt, which in turn reacts with a chelator, producing a pink-colored compound. Any intensity of the pink coloring is considered positive.

pH: Uses adual indicator stem. The red methyl and blue indicator of bromotimol are used to give a wide range of colors, covering pH levels between 5.0 and 9.0. The color scale varies from orange to greenish yellow at pH between 4.0 and 6.5, due to the color change of the red methyl indicator and from green to blue at pH between 6.5 and 9.0, due to the color change caused by the blue indicator of bromotimol.

Blood/Hemoglobin: The reaction is based on the activity of hemoglobin pseudo-peroxidase which, if present in the urine sample, catalyzes the reaction between chromogen and organic peroxide. The color scale varies from yellow, for negative samples, to bluish green, for positive samples of higher concentration.

Specific gravity: This test is based on the pKa change of the polyelectrolyte, which ionizes by releasing hydrogen ions. Low concentrations of these ions indicate low density samples, and the color action of the parameter varies from blue to greenish tones. On the other hand, high concentrations of hydrogen ions indicate high density samples, and the parameter coloration assumes yellowish tones.

Leukocytes: The reaction detects the presence of leukocyte esterases, enzymes existing in the leukocytes s. Leukocytes, if present in the urine sample, are lysed and release esterases that decompose an ester to an aromatic and acidic compound. The product of this decomposition reacts with diazonic salt, which in acidic medium produces the violet color proportional to the number of leukocytes in the urine. The color scale varies from white, when negative, to violet, when positive.

KIT PRESENTATION

REF 524050-UM: 50 Determinations

- Test strip: 50 units
- Instructions for use - available on the www.wamadiagnostica.com.br

turbidez.

6. Nitrite

Expected Values: Normalmente nitrite is not detectable in urinal.

Performance characteristics: Comparing the reagent area against a white background can aid in detecting low levels. The test is nitrite-specific and will not react with any other substance normally excreted in the urine.

Limitation of procedures: Ascorbic acid (> 25 mg/dl) may lead to a false-negative result in samples containing low concentration of nitrite in urine (<0.03 mg/dl). The negative result is not affirmative that the patient is free of bacteriuria. Pink dots or the pink s area should be interpreted as positive results. Negative results can occur when urinary tract infections are not caused by organisms that do not contain nitrate reductase, when urine has not been retained in the bladder enough (four hours or more) to reduce nitrate to nitrite to occur or when dehydrated nitrate is absent.

7. pH

Expected Values: Urine pH values usually range from 9.5 a

Performance characteristics: The test averages the pH values usually within unit 1 in the range of 5,0 – 9,0.

Limitation of procedures: Excess urine on the test strip can take acid from the buffer of the protein reagent area to the pH parameter, causing errors in reading. This phenomenon is known as “run-over”.

8. Blood/Hemoglobin

Expected Values: Hemoglobin is usually not detected in urine (10 RBC/μL; 0.03 mg/dl). When hemoglobin is detectable in the urine it may indicate kidney disease or a disorder in the urinary tract. Blood can often be found in the urine of women in menstrual period.

Performance characteristics: The test is usually able to detect 10 cells/μl with a trace.

Limitation of procedures: The test is slightly more sensitive to free hemoglobin and myoglobin than intact erythrocytes. High density or high protein can reduce blood test reactivity. Microbial peroxidase associated with urinary tract infection may present false-positive results. Concentrations of ascorbic acid (> 40 mg/dl) can cause false-negative results at low blood levels.

9. Specific gravity

Expected Values: The normal density is between 1,001 and 1,035.

Performance characteristics: The test allows the determination of density in urine between 1,000, 1,005, 1,010, 1,015, 1,020, 1,025, 1,030.

Limitation of procedures: Highly alkaline sample can lead to a decrease in the result, while highly acidic samples can significantly increase the result.

10. Leukocytes

Expected Values: Leukocytes are usually not detectable in urine.

Performance characteristics: The test is usually capable of detecting 20~25 WBC/μl with a trace.

Limitation of procedures: The test result may not always be consistent with the number of countable cells through microscopic examination. High glucose concentration, high density, high albumin level, high formaldehyde concentration, or blood presence may cause decreased test results. High concentration of oxalic acid from oxidizing agents can cause false-positive result.

QUALITY CONTROL

The test strips should be properly stored and handled before and during the test. The reaction of the reagent strips should be confirmed by testing known negative and positive samples or by multiple analytical controls containing normal and abnornal amounts of each analytical being tested.

SENSITIVITY AND SPECIFICITY

Urobilinogen: 0.1 mg/dl

Glucose: 50 - 100 mg/dl

Ketone bodies: 5 mg/dl (acetoacetic acid)

Bilirubin: 0.5 mg/dl (bilirubin)

Protein: 30 mg/dl (albumin)

Nitrite: 0.05 mg/dl (nitrite ion)

Blood: 10 RBC/ μl (0.03 mg/dl hemoglobin and erythrocyte intact).

Leukocyte: 20-25 WBC/μl (intact and lysed).

LIMITATIONS OF USE

The **Uri-Color Check** is a test to actively determine the presence of various parameters in the urine. Thus, as in all laboratory tests, definitive diagnosis or therapeutic decisions should not be based on a single result or method. Definitive clinical diagnosis should be made by the physician after analysis of clinical and laboratory data.

The effects of drugs or other metabolites on the individual tests of the test strips are not known in all cases. Therefore, it is recommended that in case of doubt, the test is repeated after withdrawal of the potential interfering agent, such as medication or vitamin supplement, etc.

The correct reading time shown on the tube label is important for the accuracy of the results, any reading outside of this will invalidate the test.

Changes in staining that appear along the margin of the test area should be ignored, a careful removal of excess urine should eliminate this phenomenon.

PRECAUTIONS AND WARNINGS

- Only for in vitro diagnostic use.
- Carefully read the instructions for use before performing the test.
- The expiration date corresponds to the last day of the month marked on the kit box label. Do not use after the expiration date.
- It should be avoided to expose the kit to high temperatures, as well as directly to the sun.
- Do not freeze the kit components as this will cause irreversible deterioration.
- Do not use the kit when any component has a visual characteristic in disagreement with the one specified in the PRODUCT FISPQ.
- Keep care used in safety when handling components. Samples should be handled as potentially infectious materials.

- Dispose of all material in accordance with local regulations. Brazil: consult Nº 222 de 28/03/18 of ANVISA.
- Use Good Laboratory Practices (PLB) in the preservation, handling and disposal of materials.

RISK MANAGEMENT

WAMA Diagnóstica, after reviewing and critically detailed analysis of all known and/or predicted hazards, concludes that all risks associated with the Uri-Color Check product were evaluated, that measures to reduce the risk were implemented and that the product does not present greater risks than the benefits obtained from its use; and that, if used by qualified and trained professionals, aware of the precautions described in the products, will perform their functions with quality, safety and efficacy.

WARRANTY TERM

WAMA Diagnóstica guarantees the exchange of this diagnostic set, if it is within the expiration date and is proven by its Technical Advisory that there were no failures in the execution, handling and conservation of this product. **WAMA Diagnóstica** and its distributors are not responsible for failures in kit performance under these conditions.

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